# Sample Requirements of MS-Services

## 1-1 /2. SDS-PAGE analysis

- (b) Sample volume:  $\leq 50 \ \mu L$
- (c) Sample concentration:  $\ge 0.2 \ \mu g/\mu L$  (Silver stain),  $\ge 1 \ \mu g/\mu L$  (Coomassie blue stain)
- (d) Platform: Hoefer SE260, SE600

## 1-3 /4. Two-dimension Electrophoresis (IEF/SDS-PAGE)

- (b) Sample volume:  $\leq 150 \ \mu L$
- (c) Sample concentration: ≥1 µg/µL (Silver stain), ≥2 µg/µL (Coomassie blue stain)
- (d) Platform: Ettan IPGphor 3 IEF system/Hoefer SE600

## 1-5. Protein solution concentration, salt removal

- (a) Sample type: Protein extracts > Purified protein solution > Body fluid
- (b) Sample amount:  $\leq 2 \text{ mL}$
- (c) Sample concentration:  $\geq 0.01 \ \mu g/\mu L$

## 1-6. Protein solution concentration, salt removal

- (a) Sample type: Protein extracts > Purified protein > Body fluid
- (b) Sample amount:  $\leq 50 \text{ mL}$
- (c) Sample concentration:  $\geq 0.01 \ \mu g/\mu L$

## 1-7. Plasma abundant proteins (Top-12) removal

- (d) Sample type: Human body fluids (serum, plasma, urine, interstitial fluid)
- (e) Sample amount:  $\geq 50 \ \mu L$
- (f) Sample concentration:  $\geq 40 \ \mu g/\mu L$

## 2-1. Purified protein: Peptide mapping

- (a) Sample type: Purified protein with amino acid sequence
- (b) Sample amount:  $\geq 20 \text{ pmol} (\sim 1 \text{ } \mu \text{g for } 50 \text{ } \text{kDa protein})$
- (c) Sample concentration:  $\geq 2 \text{ pmol}/\mu L$  (~0.1 µg/µL for 50 kDa protein)
- (d) Sample purity:  $\geq 50\%$
- (e) Platform: 1DLC, LTQ-Orbitrap MS

# 2-2. Purified protein: molecular weight determination

- (a) Sample type: Purified protein ( $\leq 70$  kDa)
- (b) Sample amount:  $\geq 100 \text{ pmol} (\sim 5 \text{ } \mu\text{g} \text{ for } 50 \text{ } \text{kDa protein})$
- (c) Sample concentration:  $\geq 10 \text{ pmol/}\mu\text{L}$  (~0.5  $\mu\text{g/}\mu\text{L}$  for 50 kDa protein)
- (d) Sample purity:  $\geq 50\%$
- (e) Platform: 1DLC, LTQ-Orbitrap MS

#### 3-1. In-gel protein identification: Basic analysis

- (a) Sample type: Gel slice (please provide gel image with staining)
- (b) Gel size:  $\leq 3 \text{ mm} \times 8 \text{ mm} \times 1.5 \text{ mm}$  (Height  $\times$  Width  $\times$  Thickness)
- (c) Platform: 1DLC, LTQ-Orbitrap MS

#### 3-2/3. In-gel protein identification: Standard/Advanced analysis

- (a) Sample type: Cell pellet, Protein extracts
- (b) Sample volume:  $\leq 50 \ \mu L$
- (c) Sample concentration:  $\ge 0.2 \ \mu g/\mu L$  (silver stain),  $\ge 1 \ \mu g/\mu L$  (Coomassie blue stain)
- (d) Platform: Hoefer SE260, 1DLC, LTQ-Orbitrap MS

#### 3-4/5/6. In-solution protein identification: Basic/Standard/Advanced analysis

- (a) Sample type: Cell pellet, protein extracts (please provide buffer constituents)
- (b) Sample volume:  $\leq 100 \ \mu L$
- (c) Sample concentration: ≥0.1 µg/µL (Basic/Standard analysis), ≥0.5 µg/µL (Advanced analysis)
- (d) Platform: 1DLC, 2DLC, LTQ-Orbitrap MS

## 4-1/2. Protein interaction: Basic/Standard analysis

- (a) Sample type: Immunoprecipitation product (beads or eluents, please provide buffer constituents for solution sample)
- (b) Sample amount: Not specified (recommended IP starting materials  $\geq$  500 µg)
- (c) Sample concentration: Not specified.
- (d) Platform: Hoefer SE260, 1DLC, LTQ-Orbitrap MS

## 4-3. Protein interaction: interacting region analysis

- (a) Sample type: Purified protein with amino acid sequence and condition of molecules interaction
- (b) Sample amount:  $\geq 200 \text{ pmol} (\sim 10 \text{ } \mu\text{g for } 50 \text{ } k\text{Da protein})$
- (c) Sample concentration:  $\geq 4.5 \text{ pmol/}\mu\text{L}$  (~0.23 µg/µL for 50 kDa protein)
- (d) Sample purity:  $\geq 90\%$
- (e) Platform: online-Pepsin-LC, Orbitrap HDX-MS

#### 5a. Quantitative proteome: Basic/Standard/Advanced analysis

- (a) Sample type: Cell pellet, Protein extracts (please provide buffer constituents)
- (b) Sample amount: ≥5 μg (Basic), ≥10 μg (Standard), ≥20 μg (Advanced) for each group
- (c) Sample concentration: ≥0.5 µg/µL (Basic), ≥1 µg/µL (Standard), ≥2 µg/µL (Advanced)
- (d) Limitations: Detergents ( $\leq 0.2\%$ ), Urea ( $\leq 1$  M); Buffer without Tris-Base, Ammonium Salt, Gly-gly, DTT, Mercaptoethanol
- (e) Platform: online-2DLC LTQ-Orbitrap MS

#### 5b. Quantitative phosphoproteome: Basic/Standard/Advanced analysis

- (a) Sample type: Cell pellet, Protein extracts (please provide buffer constituents)
- (f) Sample amount: ≥10 µg (Basic), 40 µg (Standard), 100 µg (Advanced) for each group
- (b) Sample concentration: ≥0.5 µg/µL (Basic), ≥1 µg/µL (Standard), ≥2 µg/µL (Advanced)
- (c) Limitations: Detergents ( $\leq 0.2\%$ ), Urea ( $\leq 1$  M); Buffer without Tris-Base, Ammonium Salt, Gly-gly, DTT, Mercaptoethanol
- (d) Platform: online-2DLC LTQ-Orbitrap MS

#### 6. Protein post-translational modification site analysis

- (a) Sample type: Gel slice (please provide gel image with staining), Immunoprecipitation product (beads or eluent), Protein extracts or Purified protein solution (please provide buffer constituents for solution sample)
- (b) Gel size:  $\leq 3 \text{ mm} \times 8 \text{ mm} \times 1.5 \text{ mm}$  (Height  $\times$  Width  $\times$  Thickness)
- (c) Sample amount:  $\geq 100 \ \mu g$  (Protein extracts),  $\geq 10 \ \mu g$  (Purified protein solution, purity  $\geq 50\%$ ),  $\geq 500 \ \mu g$  (recommended IP starting materials)
- (d) Sample concentration:  $\geq 0.25 \ \mu g/\mu L$
- (e) Platform: PTM-peptide enrichment, 1DLC, LTQ-Orbitrap MS

## 7. Protein quantification analysis

- (a) Sample requirement: Please contact us for further details.
- (b) Platform: Bio-Plex System

## 8-1. Bioinformatics analysis: Metabolic & Signaling pathway analysis

- (a) Sample type: Quantitative data of genomic or proteomic analysis
- (b) Platform: Please contact us for further details.

## 8-2. Bioinformatics analysis: Customized database searching

- (a) Sample type: MS raw data, Processed MS data (mzML, mzXML, mgf)
- (b) Platform: Please contact us for further details.

## 9. Rush service

(a) Service type: Available for our services with turnaround time  $\geq 2$  weeks.

# **10. Customized service**

- (a) Sample requirement: Please contact us for further details.
- (b) Platform: Please contact us for further details.